CZ – Kinesio páska

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EN – Kinesio Tape

Safety Instructions:

Before applying, ensure that the skin is dry, clean, and free of oils or creams. Do not use the tape on open wounds, irritated, or damaged skin. Follow the instructions of a physiotherapist or specialist during application to ensure the tape is applied correctly and serves its purpose. Ensure the tape is properly stretched – excessive tension may cause discomfort. The product contains hypoallergenic adhesive, but in case of irritation, redness, or itching, remove the tape immediately and wash the skin. If irritation or pain persists, seek medical advice. The tape can be worn for several days (typically 3–5 days), but regularly check for signs of skin irritation. If discomfort occurs, remove the tape. The product is water-resistant, allowing use during bathing, showering, or sauna sessions. Avoid excessive rubbing of the tape when drying to extend its lifespan. The tape is designed for single use only. After removal, it will not adhere properly to ensure functionality. Do not reuse the same tape to prevent infections or reduced effectiveness. Store the tape in a dry place at room temperature, away from direct sunlight. Keep out of reach of children. If you are unsure about the correct use of the tape or have health issues, consult a specialist. The tape does not replace medical care and should be used as a supplement to treatment.

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